



LP — 20

HYGIENE

AUSTRIA

Mouth-Nose Mask, Type-II

bacterial filtration efficiency*
99%

1. Made in Austria
2. Certified
3. 3 layers
4. High breathability due to high quality non-woven materials



1. Wash your hands using soap and water or a hand sanitizer.



2. Pull the loops of the mask back to fit them around your ears.



3. Lightly press the metal piece in your mask across your nose to fit.



4. Pull the upper and lower sides of your mask over nose and chin.



5. Avoid touching the mask while wearing.

* according to EN 14683:2019

Mask Specifications

Brand	- Hygiene Austria LP GmbH
Name	- Mouth-Nose Mask, Type-II
Type	- Non-Medical
Quantity/Units	- 50 pc. per unit
Colour	- Blue-white, upon request
Instructions and Durability	- Shelf life: 5 years - When humid, the mask is to be removed and disposed of - Non-reusable, disposable
Storage Conditions	- Store above 0°C and below 40°C in a dry place without direct sunlight
Size	- 175 x 95 mm (slight deviations possible)
Weight	- 6,21 g (slight deviations possible)
Certificates	- EN 14683:2019, Type-II - CE certified machines - TÜV tested - OEKO-TEX® Standard 100, Certificate #VN020 168438 OETI
Product specifications	- Made out of 3 layers of non-woven material, 2 outside layers spunbond nonwovens and 1 inner layer melt-blown nonwovens polypropylene - 8cm rubber-coated, extra strong and flexible nose piece

Outer Carton Specifications

Measurements Carton	- 580 x 310 x 400 mm
Quantity	- 36 units, 50 pieces each
Material	- Corrugated cardboard

Please check the authenticity of all certificates offered to you by suppliers through <https://www.eu-esf.org/covid-19/4513-covid-19-suspicious-certificates-for-ppe>
Your health will thank you.



Packaging

Packaging available in 50 piece format.

Through a light perforation, the masks can be easily pulled out from the top of the packaging unit.



* according to EN 14683:2019

Hygiene Austria LP GmbH | Palmersstrasse 6 – 8 | 2351 Wiener Neudorf | FN: 530871 v

